

Confidential Report

Our Ref: 67237







Unit 11 Westpoint Enterprise Park, Clarence Avenue,

Trafford Park, Manchester, M17 1QS, UK. Telephone: +44 (0) 161 876 4211

Email: enquiries@shirley.co.uk
Website: www.shirley.co.uk

Date: 23/11/22
Our Ref: 67237
Your Ref: 83A20038
Client: Camira Fabrics
Page: 1 of 3

Client: Camira Fabrics

Address: The Watermill

Wheatley Park

Mirfield

West Yorkshire WF14 8HE

<u>Job Title</u>: Antibacterial testing

Clients Order Ref: 83A20038

<u>Date of Receipt:</u> 18/11/2022

<u>Description of sample(s)</u>: One sample received labelled as:

1.- Manila. Batch 509403. Colour: Cambric

The sample was received in good condition.

Work Requested: AATCC 147:2011 (2016)e

This report relates only to the items tested. Enquiries concerning this report should be address to Customer Services.





Copyright © 2021 Shirley Technologies Limited. All rights reserved.



Unit 11 Westpoint Enterprise Park, Clarence Avenue, Trafford Park, Manchester, M17 1QS, UK.

> Telephone: +44 (0) 161 876 4211 Email: enquiries@shirley.co.uk

Website: www.shirley.co.uk

Date: 23/11/22
Our Ref: 67237
Your Ref: 83A20038
Client: Camira Fabrics
Page: 2 of 3

Introduction

One sample was received labelled as above, for testing to AATCC 147: 2011 (2016)e. The client requested testing the face only. Testing commenced on the 22/11/2022.

Methodology: AATCC 147:2011 (2016)e

An inoculum of *S. aureus* ATCC 6538 (or equivalent) was prepared by transferring 1 ± 0.1 ml of a 24-hour nutrient broth culture into 9 ± 0.1 ml of sterile distilled water.

Using a 4mm inoculating loop, one loopful of the diluted inoculum was transferred to the surface of a sterile nutrient agar plate by making five 60 mm long parallel streaks spaced 10 mm apart without refilling the loop.

The non-sterile test specimen was then placed transversely across the five inoculum streaks. This was performed in triplicate.

The replicates was then incubated at $37 \pm 2^{\circ}$ C for 24 hours and repeated with *K.pneumoniae* ATC 4352 (or equivalent).

The zone of inhibition around the specimen was determined according to the method stated in Section 11 in AATCC Test Method 147:2011 (2016)e.

Controls were performed ensuring that both bacteria did grow adequately.

Results

The replicate samples were examined for any zone of inhibition and if there was any growth under the sample.

The following key is used in the tables below:

W = Width of clear zone of inhibition.

S.aureus ATCC 6538

Sample: 1 Manila	Specimen	Average W (cm)	Growth under sample
Face	1	0.27	No Growth
	2	0.31	No Growth
	3	0.31	No Growth





1066



Unit 11 Westpoint Enterprise Park, Clarence Avenue, Trafford Park, Manchester, M17 1QS, UK.

> Telephone: +44 (0) 161 876 4211 Email: enquiries@shirley.co.uk Website: www.shirley.co.uk

Date: Our Ref: 23/11/22

Your Ref:

67237 83A20038

Client: Page: Camira Fabrics 3 of 3

K.pneumoniae ATCC 4352

Sample: 1 Manila	Specimen	Average W (cm)	Growth under sample
Face	1	0.35	No Growth
	2	0.33	No Growth
	3	0.39	No Growth

Comments and Conclusions

There was no contamination of the sample as this was verified with incubating under the same test conditions, without the test bacteria.

The control samples performed as expected and proved the validity of the testing.

The criterion for passing the test must be agreed upon by the interested parties.

The standard states that to constitute acceptable antibacterial activity, there must be no bacterial colonies directly under the sample in the contact area.

Where required to make a judgement to any pass/fail criteria an estimation of uncertainty of measurement has been taken into account.

Under our Policy we have used a non-binary decision rule. See our decision rules Policy (http://www.bttg.co.uk/decision-rules-policy) for further information.

Estimated uncertainty of measurement where k=2 (approximately 95% confidence level) = \pm 0.25 log units.

Reported by:

Cristina Alepuz Section Leader

Countersigned by: ...

Gareth Heywood

Copyright © 2021 Shirley Technologies Limited. All rights reserved.

Manager



